

BODE X-Wipes Safety Pack

Your versatile helper that ensures maximum compliance in high risk areas with simultaneously convenient and easy use.



Advantages at a glance

- BODE X-Wipes fleece roll in the Safety Pack with integrated filling/dispensing opening
- Flexibility known from the BODE X-Wipes system combined with a maximum of hygienic safety, as no reprocessing is necessary
- Targeted use in high-risk areas
- Good stability
- High-quality PET fleece guarantees optimal release of active ingredients
- Suitable for all liquid surface disinfectants or exclusive products for instrument reprocessing from HARTMANN

Areas of application

BODE X-Wipes Safety Pack is recommended for the disinfecting surface cleaning in all areas of health care facilities and laboratories that require maximum hygiene, cost-effectiveness, and user convenience, for example:

- on work surfaces
- on patient stretchers
- on near-patient surfaces
- on surfaces in isolation rooms
- to prevent cross contamination in critical hygiene areas
- in areas that have no water connection

BODE X-Wipes Safety Pack is especially for high risk areas, however it is also recommended for medium to low risk areas without worrying about reprocessing of dispenser.

Application


BODE X-Wipes Safety Pack can be used with all liquid surface disinfectants or exclusive products for instrument reprocessing from HARTMANN. The system is filled with 2.0 litres of disinfection solution.

When preparing use-solutions from concentrates (e.g. Mikrobac® forte, Kohrsolin® FF), choose at least the use concentration of 0.5% or higher. When preparing use-solutions, always observe the application and safety information of the respective product. Alcohol-based rapid disinfectant (e.g. Bacillol® AF, Bacillol® 30 Sensitive Foam) are ready-to-use and are used undiluted.



BODE X-Wipes Safety Pack


Filling the BODE X-Wipes Safety Pack




Put on gloves.




Unscrew lid.



Add 2.0 litres of disinfection solution.



Pull wipe through dispensing opening.



Screw lid back on Safety Pack.



Standing time:
28 days/ 7 days.

No reprocessing of the BODE X-Wipes Safety Pack necessary

Dosage

Add 2.0 litres use-solution of the needed disinfectant

Composition

Absorbent polyester fleece 60 g/m².

Compatibility

The BODE X-Wipes Safety Pack can be used with all liquid HARTMANN surface disinfectants as well as the following instrument disinfectants: Bomix[®] plus^{a)} and Bodedex[®] forte^{b)}.

In case of multiple use: Danger of contamination and mixing of products/ batches.

Standing time

28 days for surface disinfectants.
(7 days for Bodedex[®] forte).

Expert assistance

HARTMANN SCIENCE CENTER

Tel: +49 (0) 40 - 54 00 6 -111

E-Mail: science-center@hartmann.info



**Research for
infection protection**
hartmann-science-center.com

Product

BODE X-Wipes Safety Pack

Content

Dispenser with 90 wipes

Item no

on request

Please note: the availability of products in the BODE X-Wipes range may vary in different countries and regions. Contact your local distribution partner for more information. The recommendations regarding our preparations are based on scientific tests and are given in good faith. More detailed recommendations, e.g. regarding material compatibility, are possible only in separate, individual cases. Our recommendations are not binding and do not constitute a guarantee. They do not preclude a company's own testing for the intended purpose and process. In this respect we cannot accept any liability. This is in accordance with our general conditions of sale and supply.

Use disinfectants safely. Always read the label and product information before use.

(a) 0482

(b)



Helps. Cares. Protects.

PAUL HARTMANN AG
Paul-Hartmann-Straße 12
89522 Heidenheim
Germany

Phone: +49-7321-36-0
Fax: +49-7321-36-3636
info@hartmann.info
www.hartmann.info

Manufacturer: BODE Chemie GmbH, Melanchthonstraße 27, 22525 Hamburg
A company of the HARTMANN GROUP