



Bacillo® Oxy Tissues

Ready-to-use, oxygen-active disinfectant wipes in a practical jumbo-canister for disinfection of medical devices.



Characteristics

- Ready-to-use disinfection wipes
- Sporicidal effective including *C. difficile* spores
- Solution based on hydrogen peroxide
- Contains no peracetic acid, colourants, or fragrances
- Good material compatibility also suitable for sensitive surfaces and medical devices

Areas of application

Bacillo® Oxy Tissues are for cleaning disinfection of surfaces of non-invasive and invasive medical devices with enhanced efficacy. Not suitable as the end point of processing of invasive medical devices.

Microbiology

Bactericidal, yeasticidal, tuberculocidal, mycobactericidal, fungicidal, effective against *C. difficile* spores, sporicidal, limited spectrum virucidal activity, virucidal

Proven efficacy	Testing method	Condition	Exposure time
Bacteria / Fungi			
Bactericidal / Yeasticidal	EN 13727 / EN 13624 / EN 16615	clean	1 min
Tuberculocidal / Mycobactericidal	EN 14348 / EN 16615	clean	5 min
Fungicidal	EN 13624 / EN 16615	clean	15 min
Bacterial spores			
Sporicidal activity against <i>C. difficile</i>	EN 17126 / EN 17846	clean	15 min
Sporicidal	EN 17126 / EN 17846	clean	15 min
Viruses			
Virucidal against enveloped viruses (incl. HBV, HIV, HCV)	EN 14476	clean	1 min
Virucidal	EN 14476	clean	5 min

Bacillo[®] Oxy Tissues



Directions for use

Thoroughly wipe surfaces with Bacillo[®] Oxy Tissues. Ensure complete coverage so that optimum disinfection is guaranteed.

To prevent the wipes from drying out, reseal the wipe packaging immediately after use.

Dispose of used tissues.

Please refer to the manufacturer's instructions for the surfaces to be reprocessed. Do not allow any disinfectant solution to enter the interior of electrical devices

The product is not suitable for surfaces made of non-ferrous metals (like copper, brass) or marble.

Protective gloves are required. Do not use for skin cleansing.

Shelf life after opening

In opened and resealed jumbo-canister: 28 days.

Hazardous Classification



Warning

H315: Causes skin irritation.

H319: Causes serious eye irritation.

Composition

Active substances in the soaking solution:

70 mg/g Hydrogen peroxide

21 mg/g Glycolic acid

Chemical-physical data

pH value: 2.2

Flashpoint: > 100 °C

Density (20 °C): approx. 1.03 g/cm³

Refractive index: 1.34

Expert assistance

HARTMANN SCIENCE CENTER

Tel: +49 (0) 40 - 54 00 6 -111

E-Mail: science-center@hartmann.info



**Research for
infection protection**
hartmann-science-center.com

Product name

Content

Item no

Bacillo[®] Oxy Tissues

100 Tissues Flowpack
Size 20 x 20 cm

on request

Please note: That the availability of Bacillo[®] products may vary in different countries and regions. Contact your local distribution partner for more information. The recommendations regarding our preparations are based on scientific tests and are given in good faith. More detailed recommendations, e.g. regarding material compatibility, are possible only in separate, individual cases. Our recommendations are not binding and do not constitute a guarantee. They do not preclude a company's own testing for the intended purpose and process. In this respect we cannot accept any liability. This is in accordance with our general conditions of sale and supply.



Helps. Cares. Protects.

PAUL HARTMANN AG
Paul-Hartmann-Straße 12
89522 Heidenheim
Germany

Phone: +49-7321-36-0
Fax: +49-7321-36-3636
info@hartmann.info
www.hartmann.info

Manufacturer: BODE Chemie GmbH, Melanchthonstraße 27, 22525 Hamburg
A company of the HARTMANN GROUP