

HARTMANN



Mikrobac® Tissues

Alcohol-free disinfection wipes in the convenient flow pack. For the cleaning disinfection of alcohol-sensitive surfaces and surfaces of non invasive medical devices.



Advantages at a glance

- ready-to-use disinfection wipes
- easy and safe dispensing of individual wipes out of resealable flow pack
- alcohol, aldehyde, colourant and fragrance-free
- patented system solution: optimal release of active ingredients during disinfection
- also suitable for alcohol-sensitive materials
- can be used alternating with products based on other active substances (e.g. aldehydes or amines) without any problems

Intended Use

Disinfection of non-invasive medical devices

Mikrobac® Tissues

Application

Mikrobac® Tissues are suitable for the convenient cleaning disinfection of sensitive surfaces and surfaces of medical devices.

Areas of application

Mikrobac® Tissues are suitable for the uncomplicated cleaning disinfection of:

Medical equipment, such as:

- alcohol-sensitive application parts that cannot be submerged in liquids, such as ultrasonic heads of probes for abdominal examinations.

All kinds of alcohol-sensitive medical inventory and surfaces, such as:

- surfaces in hospitals, medical and dental practices, ambulances
- surfaces in medical laboratories, residential and nursing homes
- small objects and surfaces such as toilet seats, door handles, bed frames, and office desks

Directions for use

Carefully wipe surfaces with Mikrobac® Tissues. Ensure complete coverage so that optimum disinfection is guaranteed. If necessary, wipe off surface with a single use cloth after expiry of exposure time (e.g. if people are likely to come into direct contact with the surface). Dispose of tissue appropriately after use.

To remove disinfectant residue of medical devices, wipe surfaces with a cloth soaked in water (at least drinking water quality) after the exposure time.

We recommend wearing suitable gloves while using this product. Not to be used for skin cleansing.

When switching from one product to another, an intermediate cleaning has to be carried out.

Please follow the manufacturer's instructions. Do not allow disinfection solution to get inside of electrical devices. Not to be used for the final disinfection of semicritical medical devices (e.g. ultrasonic heads of probes for vaginal examinations).

In order to prevent the wipes from drying out close lid when not in use.

Microbiology

- bactericidal
- yeasticidal
- virucidal against enveloped viruses (incl. HBV, HIV, HCV)
- MNV
- rota- and polyomavirus

Composition

Active ingredients of the solution:

Benzyl-C12-18-alkyldimethyl-ammonium chloride 4 mg/g;
didecyldimethylammonium chloride 4 mg/g




Mikrobac® Tissues

Directions for use and Dosage

Mikrobac® Tissues are ready-to-use.

Efficacy	Testing method	Condition	Exposure time
Bacteria and Fungi			
Bactericidal / Yeasticidal	EN 16615 / EN 13727 / EN 13624	dirty	30 sec
	VAH*	dirty	1 min
Viruses			
Virucidal activity against enveloped viruses (incl. HBV, HIV, HCV)	EN14476	clean + dirty	15 sec
	DVV		30 sec
Polyomavirus	DVV		1 min
Rotavirus	DVV		30 sec
Norovirus	EN 14476	clean + dirty	4 h
Food			
Bactericidal / Yeasticidal food processing area	EN 13697 / EN 1276 / EN 1650	clean (4 °C, 10 °C and 20 °C)	1 min
		dirty (4 °C, 10 °C and 20 °C)	5 min



* mechanical

Recommendations are highlighted

Chemical-physical data

All data referring to the wipes:

Dimensions:

Flow pack with 80 tissues: 180 x 200 mm

Flow pack with 40 XXL tissues: 250 x 380 mm

Reach:

Tissue 180 x 200 mm: approx. 1.98 m²

Tissue 250 x 380 mm: approx. 2.7 m²

The reach of Mikrobac® Tissues not only depends on their size, but also on further factors, for example, the ambient temperature and structure of the surface to be disinfected. Always make sure to completely cover the surface area.

Material Compatibility

Metals:

Stainless steel (V2A, V4A), not compatible with aluminium, copper, brass

Plastics:

Polyamide (PA), polyethylene (PE), polypropylene (PP), acrylonitrile butadiene styrene copolymer (ABS), polystyrene (PS), polyvinyl chloride (PVC), rubber, polycarbonate (Macrolon), PMMA (plexi glass), polytetrafluoroethylene (PTFE), polysulfone (PSU), polyurethane (PUR), latex, silicone

Mikrobac® Tissues

Product name	Content	Item no
Mikrobac® Tissues	80 tissues flow pack size 180 x 200 mm	on request
Mikrobac® Tissues in XXL	40 tissues flow pack size 250 x 380 mm	on request
Flowpack Wallholder	1 wallholder	on request
Flowpack-holder flex	1 holder	on request



Please note: That the availability of Mikrobac® products may vary in different countries and regions. Contact your local distribution partner for more information. The recommendations regarding our preparations are based on scientific tests and are given in good faith. More detailed recommendations, e.g. regarding material compatibility, are possible only in separate, individual cases. Our recommendations are not binding and do not constitute a guarantee. They do not preclude a company's own testing for the intended purpose and process. In this respect we cannot accept any liability. This is in accordance with our general conditions of sale and supply.

*Use disinfectants safely. Always read the label and product information before use.
Please also read instructions for use and attached instructions.*



Helps. Cares. Protects.

Manufacturer:
BODE Chemie GmbH
A company of the HARTMANN GROUP
Tel: +49 40 5 40 06-0 · Fax: -200
info@bode-chemie.de
www.bode-chemie.com