



Mikrobac® Tissues / -XXL

Alcohol-free cleaning and disinfecting wipes for alcohol-sensitive surfaces in a resealable flow pack.



CE 0482

Characteristics

- Cleaning and disinfecting in one step
- Rapidly effective
- Colourant-, perfume-, alcohol- and aldehyde-free
- Good material compatibility
- Can be used alternating with products based on other active substances (e.g. aldehydes or amines) without any problems

Areas of application

Mikrobac® Tissues are QAV-based ready-to-use wipes for the cleaning and disinfection of surfaces of non-invasive medical devices. Suitable for alcohol-sensitive surfaces and e.g. for the cleaning and disinfection of non-critical ultrasound heads.

For larger surfaces, the use of Mikrobac® Tissues in XXL format is recommended:

- e.g. patient couches, operating/x-ray tables, wheelchairs/aids, surfaces in ambulances etc.

Microbiology

Bactericidal, yeasticidal, limited spectrum virucidal activity (incl. HBV, HIV, HCV), norovirus, polyoma and rotavirus

Efficacy	Testing method	Condition	Exposure time
Bacteria / Fungi			
Bactericidal / Yeasticidal	EN 16615/ EN 13727/ EN 13624	dirty	30 sec
Viruses			
Virucidal activity against enveloped viruses (incl. HBV, HIV, HCV)	EN 14476	dirty	15 sec
Norovirus	EN 14476	dirty	4 h
Food / industry			
Bactericidal / Yeasticidal	EN 13697/ EN 1276/ EN 1650	dirty (4 °C, 10 °C und 20 °C)	5 min

Mikrobac® Tissues / -XXL



Directions for use

Wear gloves when handling the product. Wipe all surfaces. Ensure complete coverage. To remove disinfectant residues after the exposure time, wipe the surfaces, especially surfaces that come into contact with food, with water (at least drinking water quality). Reseal the package directly after use. Intermediate cleaning must be carried out when changing products. Disinfectable area (per tissue)
- 80 tissues (18 cm x 20 cm): 1.5 m²
- 40 tissues XXL (25 cm x 38 cm): 2.5 m²
Do not treat food and feed directly. Treat only smooth (non-absorbent / non-porous) materials. Notes: Ensure that disinfectant does not penetrate electronic equipment. Follow the manufacturer's instructions for reprocessing. Not to be used for skin cleansing.

Stability

In open and closed flowpack: 3 months.

Hazardous Classification

H412: Harmful to aquatic life with long lasting effects.

Composition

Benzyl-C12-18-alkyldimethylammonium chloride 4mg/g;
didecyldimethylammonium chloride 4mg/g

Chemical-physical data

The information refers to the

Mikrobac® Tissues and Mikrobac® Tissues XXL soaking solution:

Density (20 °C) approx. 1 g /cm³

PH Value (20 °C) approx. 8

Expert assistance

HARTMANN SCIENCE CENTER

Tel: +49 (0) 40 - 54 00 6 -111

E-Mail: science-center@hartmann.info



**Research for
infection protection**
hartmann-science-center.com

Product name	Content	Item no
Mikrobac® Tissues	80 tissues flow pack size 18 x 20 cm	on request
Mikrobac® Tissues in XXL	40 tissues flow pack size 25 x 38 cm	on request
Flowpack wallholder	1 wallholder	on request
Flowpack-holder flex	1 holder	on request

Please note: That the availability of Mikrobac® products may vary in different countries and regions. Contact your local distribution partner for more information.

The recommendations regarding our preparations are based on scientific tests and are given in good faith. More detailed recommendations, e.g. regarding material compatibility, are possible only in separate, individual cases. Our recommendations are not binding and do not constitute a guarantee. They do not preclude a company's own testing for the intended purpose and process. In this respect we cannot accept any liability. This is in accordance with our general conditions of sale and supply.

Use surface disinfectants safely. Always read the label and product information before use.



Helps. Cares. Protects.

PAUL HARTMANN AG
Paul-Hartmann-Straße 12
89522 Heidenheim
Germany

Phone: +49-7321-36-0
Fax: +49-7321-36-3636
info@hartmann.info
www.hartmann.info

Manufacturer: BODE Chemie GmbH, Melanchthonstraße 27, 22525 Hamburg
A company of the HARTMANN GROUP