

HARTMANN



Sterillium® pure

The first allergy-friendly liquid hand disinfectant for sensitive skin - certified by ECARF. Featuring the proven Sterillium® skincare complex and free of colorants and fragrances.



Image may vary according to local requirements

Application

Sterillium® pure is used as ready-to-use alcohol-based hand disinfection product – independently of water and wash-basin – to prevent infection in all areas of healthcare and industry where hygiene is important as well as in home dialysis and when travelling. Sterillium® pure contains no colourants and fragrances. Thus, it is particularly well-suited for users suffering from allergies, hypersensitivities, atopic skin and asthmatics which is proven by the ECARF certificate.

Directions for use

Hygienic hand disinfection: Apply sufficient amount (min. 3 ml) of undiluted hand disinfectant to the cupped hand so that all areas of the hands can be wetted thoroughly and completely. Carefully rub in hand disinfectant for 30 seconds – ensure to cover all areas of skin.

Characteristics

- comprehensive efficacy, including limited spectrum virucidal activity
- inactivates Norovirus within the hygienic hand disinfection
- ECARF-certified and evidentially tested: The liquid anti-allergic hand disinfectant for the sensitive skin, atopic skin and asthmatics
- suitable for sensitive skin due to absence of fragrances & colorants
- increases skin hydration, preserves skin elasticity and the natural smoothness of the skin
- suitable also for the food industry (HACCP evaluated)



Proven efficacy	Testing method	Exposure time
Bacteria / Fungi		
Hygienic hand disinfection	VAH / EN 1500	30 sec
Surgical hand disinfection	VAH / EN 12791	1.5 min
Bactericidal	EN 13727	15 sec
Yeasticidal	EN 13624	15 sec
Tuberculocidal	EN 14348	15 sec
Mycobactericidal	EN 14348	15 sec
Viruses		
Virucidal activity against enveloped viruses (incl. HBV, HIV, HCV)	DVV / EN 14476	15 sec
Limited spectrum virucidal activity	EN 14476	1 min
	prEN 17430	30 sec
Norovirus	EN 14476	15 sec
Rotavirus	EN 14476	15 sec

Surgical hand disinfection: Thoroughly rub Sterillium® pure into hands and forearms up to the elbows. When

treating hands, pay particular attention to fingertips, palms and thumbs. Hands and forearms are to be kept wet during the entire application time of 1.5 minutes.

Sterillium® pure



Stability

In opened, and tightly closed container, with mounted single-use pump or when using the Eurodispenser 1 plus / 1 plus touchless: 12 months.

Expert assistance

HARTMANN SCIENCE CENTER

Tel: +49 (0) 40 - 54 00 6 -111

E-Mail: science-center@hartmann.info



Research for
infection protection
hartmann-science-center.com

Composition

Active agents per 100 g

Propan-2-ol 45.0 g, Propan-1-ol 30.0 g

Microbiology

Bactericidal, yeasticidal, tuberculocidal, mycobactericidal, virucidal activity against enveloped viruses, limited spectrum virucidal activity.

Chemical-physical data

Appearance transparent, colorless

Density (at 20°C) approx. 0.85 g/cm³

Flashpoint 23°C
(acc. to DIN 51755)

Product	Content	Item no.
Sterillium® pure	100 ml bottle	on request
	500 ml bottle	on request
	1000 ml bottle	on request
	1000 ml CleanSafe up-side-down bottle with pre-mounted dosing pump	on request
	5 L canister	on request

Please note: that the availability of products in the Sterillium® range may vary in different countries and regions. Contact your local distribution partner for more information. The recommendations regarding our preparations are based on scientific tests and are given in good faith. More detailed recommendations, e.g. regarding material compatibility, are possible only in separate, individual cases. Our recommendations are not binding and do not constitute a guarantee. They do not preclude a company's own testing for the intended purpose and process. In this respect we cannot accept any liability. This is in accordance with our general conditions of sale and supply.

Use disinfectants safely. Always read the label and product information before use.

Please amend in accordance with local requirements (e.g. law of advertising, product status, CLP labelling)



Helps. Cares. Protects.

Manufacturer:

BODE Chemie GmbH

A company of the HARTMANN GROUP

Tel: +49 40 5 40 06-0 · Fax: -200

info@bode-chemie.de

www.bode-chemie.com